



Chromab. lid SPE with gasket f. 24 col.

Macherey Nagel

CHROMABOND lid (PP) SPE with gasket, for 24 columns (inclusive 24 Luer fittings for lid and 24 valves (plastic)) suitable for REF 730151N



Code	Description	Packaging
MNL730176N	Chromab. lid SPE with gasket f. 24 col.	1 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

