



**Colonna Chromab. SiOH, 3 mL, 200 mg**

**Macherey Nagel**

CHROMABOND column SiOH volume: 3 mL, content of sorbent: 200 mg  
material: glass, with glass fibre filter element pack of 50



| Code       | Description                         | Packaging |
|------------|-------------------------------------|-----------|
| MNL730214G | Chromab. columns SiOH, 3 mL, 200 mg | 50 pz     |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

