



**Chromab. tubing adaptor f. 15, 45, 70 mL**

**Macherey Nagel**

CHROMABOND tubing adaptor for application of large sample volumes to 15, 45 and 70 mL PP columns incl. PTFE tube (approx. 1 m) with weight pack of 4



Code	Description	Packaging
MNL730386	Chromab. tubing adaptor f. 15, 45, 70 mL	4 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

