



Adattatore per tubi Chromab. 3, 6 mL glass

Macherey Nagel

CHROMABOND tubing adaptor for application of large sample volumes to glass columns (3 mL, 6 mL) incl. PTFE tube (approx. 1 m) with weight pack of 4



| Code | Description | Packaging |
|-----------|--|-----------|
| MNL730387 | Chromab. tubing adaptor f. 3, 6 mL glass | 4 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

