



Cartuccia filtrante CHROMAFIL PTFE, 3 mL, 0.45 μ m

Macherey Nagel

CHROMAFIL filtration cartridges PTFE volume: 3 mL membrane material:
polytetrafluoroethylene (PTFE) pore dia.: 0.45 μ m pack of 100



Code	Description	Packaging
MNL730570345	CHROMAFIL filtr. cart. PTFE, 3 mL,0.45 µm	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

