

**Colonna Chromab. HR-XAW (45 μ m),3mL,200mg****Macherey Nagel**

CHROMABOND columns HR-XAW (45 μ m, fine grained) weak mixed-mode polymer-based anion exchanger particle size: 45 μ m volume: 3 mL, content of sorbent: 200 mg material: PP, with PE filter elements pack of 30



Code	Description	Packaging
MNL730748P45	Chromab. columns HR-XAW (45 µm),3mL,200mg	30 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

