

**Chromab. starter kit Flash****Macherey Nagel**

CHROMABOND Flash starter kit consist of: 1 Tube, material: PTFE, OD: 1/8", ID: 1,5 mm, length: 3 m 5 Nut, material: PP, connection: 1/4"-28 5 Ferrules, material: ETFE, connection: 1/8" 5 Union, material: Nylon, connection: 1/4"-28 2 Female Luer lock, material: PP, connection: 1/4"-28 1 Male Luer lock, material: PP, connection: 1/4"-28 1 Male Luer, material: PP, connection: 1/4"-28



Code	Description	Packaging
MNL730798	Chromab. starter kit Flash	1 kit



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

