



Colonna Chromab. HR-XA, 3 mL, 200 mg

Macherey Nagel

CHROMABOND columns HR-XA strong mixed-mode polymer-based anion exchanger volume: 3 mL, content of sorbent: 200 mg material: PP, with PE filter elements pack of 30



Code	Description	Packaging
MNL730951	Chromab. columns HR-XA, 3 mL, 200 mg	30 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

