



Chromafix HR-XAW (S)

Macherey Nagel

Solid Phase Extraction Cartridges CHROMAFIX HR-XAW (S) - weak mixed-mode polymer-based anion exchanger - minimum adsorbent weight: 50 mg pack of 50



Code	Description	Packaging
MNL731771	Chromafix HR-XAW (S)	50 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

