



Chromafix HR-XCW (L)

Macherey Nagel

Solid Phase Extraction Cartridges CHROMAFIX HR-XCW (L) weak mixed-mode polymer-based cation exchanger minimum adsorbent weight: 450 mg pack of 50



| Code | Description | Packaging |
|-----------|----------------------|-----------|
| MNL731776 | Chromafix HR-XCW (L) | 50 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

