



Chromafix PS-Mix (S) (PS-H+/PS-OH-)

Macherey Nagel

Solid Phase Extraction Cartridges CHROMAFIX PS-Mix (PS-H+ und PS-OH-) (S) minimum adsorbent weight: 70 mg pack of 50



| Code | Description | Packaging |
|-----------|-------------------------------------|-----------|
| MNL731909 | Chromafix PS-Mix (S) (PS-H+/PS-OH-) | 50 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

