



Chromab. Flash RS 200 SiOH 15-40 μm

Macherey Nagel

CHROMABOND columns Flash RS 200 filled with: SiOH, 15-40 μm filling: 200 g
material: PP, with PE filter elements pack of 6



Code	Description	Packaging
MNL732706	Chromab. Flash RS 200 SiOH 15-40 µm	6 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

