



Chromab. Flash RS 1600 SiOH 15-40 μm

Macherey Nagel

CHROMABOND columns Flash RS 1600 filled with: SiOH, 15-40 μm filling: 1600 g material: PP, with PE filter elements pack of 2



Code	Description	Packaging
MNL732709	Chromab. Flash RS 1600 SiOH 15-40 µm	2 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

