



**Chromab. Flash RS 1600 SiOH 15-40  $\mu\text{m}$**

**Macherey Nagel**

CHROMABOND columns Flash RS 1600 filled with: SiOH, 15-40  $\mu\text{m}$  filling: 1600 g material: PP, with PE filter elements pack of 2



Code	Description	Packaging
MNL732709	Chromab. Flash RS 1600 SiOH 15-40 µm	2 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

