



Chromab. Flash RS 800 C18 ec, 15-40 μm

Macherey Nagel

CHROMABOND columns Flash RS 800 filled with: C18 ec, 15-40 μm filling: 880 g material: PP, with PE filter elements pack of 1



Code	Description	Packaging
MNL732718	Chromab. Flash RS 800 C18 ec, 15-40 µm	1 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

