



Chromab. Flash RS 80 SiOH

Macherey Nagel

CHROMABOND columns Flash RS 80 filled with: SiOH filling: 80 g Material: PP,
with PE-frits pack of 12



| Code | Description | Packaging |
|-----------|---------------------------|-----------|
| MNL732804 | Chromab. Flash RS 80 SiOH | 12 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

