



**Chromab. Flash RS 800 SiOH**

**Macherey Nagel**

CHROMABOND columns Flash RS 800 filled with: SiOH filling: 800 g Material: PP, with PE-frits  
pack of 2



Code	Description	Packaging
MNL732808	Chromab. Flash RS 800 SiOH	2 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

