



**Chromab. Flash RS 3000 SiOH**

**Macherey Nagel**

CHROMABOND columns Flash RS 3000 filled with: SiOH filling: 3000 g  
Material: PP, with PE filter elements pack of 1



| Code      | Description                 | Packaging |
|-----------|-----------------------------|-----------|
| MNL732850 | Chromab. Flash RS 3000 SiOH | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

