



Chromab. Flash BT 4 SiOH

Macherey Nagel

CHROMABOND columns Flash BT 4 filled with: SiOH filling: 4 g material: PP, with PE-frits outlet:
male Luer Lock pack of 20



Code	Description	Packaging
MNL732960	Chromab. Flash BT 4 SiOH	20 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

