



Chromab. Flash BT 25 SiOH

Macherey Nagel

CHROMABOND columns Flash BT 25 filled with: SiOH filling: 25 g material: PP, with PE-frits
outlet: male Luer Lock pack of 15



| Code | Description | Packaging |
|-----------|---------------------------|-----------|
| MNL732962 | Chromab. Flash BT 25 SiOH | 15 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

