

**Chromab. Flash BT 40 SiOH****Macherey Nagel**

CHROMABOND columns Flash BT 40 filled with: SiOH filling: 40 g material: PP,  
with PE-frits outlet: male Luer Lock pack of 15



Code	Description	Packaging
MNL732963	Chromab. Flash BT 40 SiOH	15 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

