

**Chromab. Flash BT 80 SiOH****Macherey Nagel**

CHROMABOND columns Flash BT 80 filled with: SiOH filling: 80 g material: PP,
with PE filter elements outlet: male Luer Lock pack of 12



Code	Description	Packaging
MNL732964	Chromab. Flash BT 80 SiOH	12 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

