



**Chromab. Flash BT 120 SiOH**

**Macherey Nagel**

CHROMABOND columns Flash BT 120 filled with: SiOH filling: 120 g material:  
PP, with PE filter elements outlet: male Luer Lock pack of 10



Code	Description	Packaging
MNL732965	Chromab. Flash BT 120 SiOH	10 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

