



Chromab. Flash BT 120 SiOH

Macherey Nagel

CHROMABOND columns Flash BT 120 filled with: SiOH filling: 120 g material: PP, with PE filter
elements outlet: male Luer Lock pack of 10



Code	Description	Packaging
MNL732965	Chromab. Flash BT 120 SiOH	10 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

