



Chromab. Flash DL 4, colonna vuota

Macherey Nagel

CHROMABOND Flash DL 4 Flash Direct Load cartridges - empty columns for solid injection
volumen: 8 mL (for app. 4 g silica or 3 g kieselguhr) material: PP incl. 3 PE filter elements (one
PE filter element inserted, two aside) inlet: female Luer Lock outlet: male Luer Lock pack of 50



Code	Description	Packaging
MNL732980	Chromab. Flash DL 4, empty columns	50 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

