

**Chromab. Flash DL 25, colonna vuota****Macherey Nagel**

CHROMABOND Flash DL 25 Flash Direct Load cartridges - empty columns for solid injection
volumen: 45 mL (for app. 25 g silica or 15 g kieselguhr) material: PP incl. 3 PE filter elements
(one PE filter element inserted, two aside) inlet: female Luer Lock outlet: male Luer Lock pack of
50



Code	Description	Packaging
MNL732982	Chromab. Flash DL 25, empty columns	50 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

