



Chromab. Flash DL 40, PE filter elements

Macherey Nagel

CHROMABOND PE filter elements for CHROMABOND Flash Direct Load cartridges DL 40
material: PE pack of 250



Code	Description	Packaging
MNL732983FE	Chromab. Flash DL 40, PE filter elements	250 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

