

**Chromab. Flash DL 80, empty column****Macherey Nagel**

CHROMABOND Flash DL 80 Flash Direct Load cartridges - empty columns for solid injection volumen: 160 mL (for app. 80 g silica or 60 g kieselguhr) material: PP incl. 3 PE filter elements (one PE filter element inserted, two aside) inlet: female Luer Lock outlet: male Luer Lock pack of 20



Code	Description	Packaging
MNL732984	Chromab. Flash DL 80, empty column	20 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

