

**Chromab. Flash DL 80, empty column****Macherey Nagel**

CHROMABOND Flash DL 80 Flash Direct Load cartridges - empty columns for solid injection
volumen: 160 mL (for app. 80 g silica or 60 g kieselguhr) material: PP incl. 3 PE filter elements
(one PE filter element inserted, two aside) inlet: female Luer Lock outlet: male Luer Lock pack of
20



Code	Description	Packaging
MNL732984	Chromab. Flash DL 80, empty column	20 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

