



Chromab. Flash DL 120, colonna vuota

Macherey Nagel

CHROMABOND Flash DL 120 Flash Direct Load cartridges - empty columns for solid injection volumen: 220 mL (for app. 120 g silica or 80 g kieselguhr)
material: PP incl. 3 PE filter elements (one PE filter element inserted, two aside)
inlet: female Luer Lock outlet: male Luer Lock pack of 20



Code	Description	Packaging
MNL732985	Chromab. Flash DL 120, empty columns	20 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

