

**Chromab. Flash DL 200, empty columns****Macherey Nagel**

CHROMABOND Flash DL 200 Flash Direct Load cartridges - empty columns for solid injection volumen: 410 mL (for app. 200 g silica or 150 g kieselguhr)
material: PP incl. 3 PE filter elements (one PE filter element inserted, two aside)
inlet: female Luer Lock outlet: male Luer Lock pack of 10



Code	Description	Packaging
MNL732986	Chromab. Flash DL 200, empty columns	10 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

