



VP 150/40 NUCLEODUR C18 HTec, 5 μ m

Macherey Nagel

VarioPrep HPLC column VP 150/40 NUCLEODUR C18 HTec, 5 μ m *** with adjustable Endfitting

*** length: 150 mm, ID: 40 mm pack of 1



Code	Description	Packaging
MNL762555400	VP 150/40 NUCLEODUR C18 HTec, 5 µm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

