



Piastre TLC precoattate RP-2 UV254 (gel di nsilice silanizzato) spessore: 0.25 mm

Macherey Nagel

size: 10x20 cm pack of 50



Code	Description	Packaging
MNL811081	TLC precoated plates RP-2 UV254 (silanized silica gel) thickness: 0.25 mm	50 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

