



Piastra TLC precoatata RP-2 UV254 (gel di silice silanizzato) spessore: 0.25 mm

Macherey Nagel

size: 20x20 cm pack of 25



Code	Description	Packaging
MNL811082	TLC precoated plates RP-2 UV254 (silanized silica gel) thickness: 0.25 mm	25 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

