

**Beuta in vetro boro 3.3, collo stretto****BRAND**

Beuta in vetro boro 3.3, collo stretto - secondo DIN ISO 1773 - con bordo salvagoccia e graduazione.



Code	Description	Packaging
LLG04689213	beuta in vetro - collo stretto - cap. 100 mL - dim. 64x105 h mm	10 pz
LLG04689214	beuta in vetro - collo stretto - cap. 300 mL - dim. 87x156 h mm	10 pz
LLG04689215	beuta in vetro - collo stretto - cap. 500 mL - dim. 105x180 h mm	10 pz
LLG04689216	beuta in vetro - collo stretto - cap. 1000 mL - dim. 131x220 h mm	10 pz
LLG04689217	beuta in vetro - collo stretto - cap. 2000 mL - dim. 166x280 h mm	10 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

