


**Tubo in silicone****BRAND**

Tubo in silicone - traslucido - autoclavabile (121 °C) - bassa resistenza chimica agli acidi concentrati.



Code	Description	Packaging
LLG04688874	tubo in silicone - spess. 1,5 mm - diam. int. 6 - diam. est. 9 mm, ATTENZIONE! ORDINE MINIMO 25 CONF. O MULTIPLI DI 25 CONF.	25 m



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

