


**Vetri orologio in vetro soda-calcico****BRAND**

Vetri da orologio - in vetro soda-calcico - in accordo DIN 12341 - bordi levigati - temperati per ridurre il rischio di rottura.



Code	Description	Packaging
LLG04688919	vetro d'orologio - diam. 60 mm	10 pz
LLG04688920	vetro d'orologio - diam. 80 mm	10 pz
LLG04688921	vetro d'orologio - diam. 100 mm	10 pz
LLG04688922	vetro d'orologio - diam. 120 mm	10 pz
LLG04688923	vetro d'orologio - diam. 125 mm	10 pz
LLG04689005	vetro d'orologio - diam. 40 mm	10 pz
LLG04689006	vetro d'orologio - diam. 50 mm	10 pz
LLG04689007	vetro d'orologio - diam. 150 mm	10 pz
LLG04689008	vetro d'orologio - diam. 200 mm	10 pz
LLG04689021	vetro d'orologio - diam. 90 mm	10 pz
LLG04689038	vetro d'orologio - diam. 70 mm	10 pz

