

Technical Specification

PRODUCT : Quinoline RE - Pure

CODE : 333700

MOLECULAR WEIGHT :129,16

CAS N° : 91-22-5

FORMULA : C₉H₇N

TEST	U.M.	SPECIFICATION
Description	-	Yellow -brown clear liquid
Identification	-	Positive
Density at 20°C	-	1.088 - 1.100
Refractive index at 20°C	-	1.6218 - 1.6318
Boiling point	°C	235.6 - 238.6
Assay (GC)	%	>= 95

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in