

Technical Specification

PRODUCT : Xylene, mix of isomers RS RSE For electronic use

CODE : 492350 **MOLECULAR WEIGHT :106,17**

FORMULA : C8H10

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Colour	APHA	<= 10
Identification (I.R.)	-	Positive
Readily carbonizable substances	-	Conform
Density at 20°C	-	0.864 - 0.870
Resistivity	MOhm.cm	>= 1
Water (K.F.)	ppm	<= 100
Residue on evaporation	ppm	<= 5
Acidity (HCl)	ppm	<= 5
Alcalinity (NH3)	ppm	<= 1
Benzene	ppm	<= 100
Chloride	ppm	<= 3
Phosphate	ppm	<= 1
Heavy metals (as Pb)	ppm	<= 0.1
Toluene	ppm	<= 5000
Total sulphur	ppm	<= 3
Ag	ppm	<= 0.02
Al	ppm	<= 0.05
As	ppm	<= 0.01
Au	ppm	<= 0.05
B	ppm	<= 0.01
Ba	ppm	<=0.1
Be	ppm	<= 0.02
Bi	ppm	<= 0.02
Ca	ppm	<= 0.2
Cd	ppm	<= 0.005
Co	ppm	<= 0.01
Cr	ppm	<= 0.01
Cu	ppm	<= 0.01
Fe	ppm	<= 0.05
Ga	ppm	<= 0.02
In	ppm	<= 0.02
K	ppm	<= 0.1
Li	ppm	<= 0.02
Mg	ppm	<= 0.1
Mn	ppm	<= 0.01
Mo	ppm	<= 0.05
Na	ppm	<= 0.1
Ni	ppm	<= 0.01
Pb	ppm	<= 0.01
Pt	ppm	<= 0.05
Sb	ppm	<= 0.01
Sn	ppm	<= 0.02
Sr	ppm	<= 0.02
Ti	ppm	<= 0.05
Tl	ppm	<= 0.05
V	ppm	<= 0.05

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.

CARLO ERBA Reagents w ould not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in

Technical Specification

TEST	U.M.	SPECIFICATION
Zn	ppm	<= 0.01
Zr	ppm	<= 0.05

General description :
Boiling point : 140 °C

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in