

Technical Specification

PRODUCT : DIISOPROPYL ETHER RE - Pure - for synthesis W

CODE : P04339 **MOLECULAR WEIGHT :102.17**

CAS N° : 108-20-3

FORMULA : [(CH3)2CH]2O

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification (I.R.)	-	Conform
Refractive index at 20°C	-	1.366 - 1.370
Free acid (as CH3COOH)	% (w /w)	<= 0.005
Peroxide (as H2O2)	ppm	<= 5
Water (K.F.)	% (w /V)	<= 0.1
Non volatile residue (w ithout stab.)	ppm	<= 30
Stabilizer (ionol)	ppm	3 - 100
Assay (GC)	%	>= 98.5
2-Propanol	% (w /V)	<= 0.1

To be used w ithin 12 months from the date of packaging

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents w ould not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in