

Technical Specification

PRODUCT : PIPERIDINE RE - Pure

CODE : P06602

MOLECULAR WEIGHT :85.15

CAS N° : 110-89-4

FORMULA : NH(CH₂)₄CH₂

TEST	U.M.	SPECIFICATION
Refractive index at 20°C	-	1.450 - 1.454
Identification (I.R.)	-	Conform
Density d20/4	-	0.856 - 0.866
Water (K.F.)	ppm	<= 5000
Assay (GC)	%	>= 99

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in