

Technical Specification

PRODUCT : PYRIDINE RE - Pure - for synthesis

CODE : P06702 **MOLECULAR WEIGHT :79.10**

CAS N° : 110-86-1

FORMULA : C5H5N

TEST	U.M.	SPECIFICATION
Refractive index at 20°C	-	1.508 - 1.512
Water (K.F.)	ppm	<= 1000
Residue on evaporation	ppm	<= 50
Colour	APHA	<= 20
Assay (GC)	%	>= 99.8

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in