

Technical Specification

PRODUCT : n-Octane PS Pure for synthesis

CODE : P08430 **MOLECULAR WEIGHT :114,23**

CAS N° : 111-65-9

FORMULA : CH3(CH2)6CH3

TEST	U.M.	SPECIFICATION
Refractive index at 20°C	-	1.395 - 1.399
Colour	APHA	<= 10
Assay (GC)	%	>= 99

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in