

Technical Specification

PRODUCT : MOLECULAR SIEVES 3 Å RS - Beads
CODE : P18100

TEST	U.M.	SPECIFICATION
------	------	---------------

Total volatile (950°C) ~ 1.5
H2O-Adsorption capacity (10% r.h.,25°C) ~ 17.5 %
H2O-Adsorption capacity (80% r.h.,25°C) ~ 21.5 %
Crush strength ~ 28 N
Particules size : 1.6-3.15 mm

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in