



CARLO ERBA Reagents, well known manufacturer of chemicals for more than 150 years with 2 production sites located in France, is at the service of its customers to meet their expectations by respecting their requirements and constraints.

CARLO ERBA Reagents has become for several years **THE specialist of the Tailor made service** for Industry and Laboratory. We are tailoring our production tool and quality know how to your specific needs.

PRODUCT

PACKAGING

PRODUCTION

Service

Flexibility

**TAILOR MADE**

Versatility

Expertise

QUALITY


Each customer has issues that may require the implementation of dedicated solutions. Indeed, request are numerous and unique. Our expertise in various steps of the product value chain allows us to give you an answer to your specific needs.

The proximity and experience of our Quality, Production and Logistics teams enable us to work with you to find the most suitable solution and make CARLO ERBA Reagents your **partner in choice**.

## SPECIFIC PURITY

According to your requirements CARLO ERBA Reagents can work on the purity of its regular products :

- Batch selection
- Distillation
- Purification
- New source
- Impurities identification



## OUT-SOURCING SERVICES



## SAFETY STOCKS

CARLO ERBA Reagents can keep a safety stock for you to palliate any unfortunate events and to avoid the storage of large quantity of chemicals in your facility.



## CUSTOM MIXTURES

CARLO ERBA Reagents can produce for you all kind of solutions or mixtures for various applications from 2.5L to several tons. This service allows you to save time, to reduce risks associated with the preparation of these solutions, and to receive a record label and safety data conform to current legislation.

- Specific solution mixtures (organic or aqueous)
- Buffered, acidified phases
- Analytical preparative and mobile phases
- Reagents for biotechnology
- Titrated solutions
- Reagents for on-line analysers like TOCmetres or silicemetres
- Standards for AAS, ICP, ICP-MS mono or multi-elements
- .....

## PURIFICATION

CARLO ERBA Reagents provides expertise with over 150 years of experience in the purification of your chemical products, custom recycling of your high purity solvents, grinding for powder.

- An initial technical study to validate project feasibility and establish your quotation is done.
- A pre-production trial is then carried out.
- Your products are received in drums, reusable stainless steel drums or disposable containers.
- They are purified or recycled according to your specifications and dispatched in the packaging of your choice, accompanied by their certificate of analysis.






# QUALITY

## QUALITY ASSURANCE


Our quality assurance system is based on risk management of technical systems and production sites, qualification of facilities adapted to our business, process control, cleaning procedures and analytical methods.

Independent of the production, the quality department manages the documentation, the control of the recordings, the traceability, the internal audits, the change control, the periodic monitoring of the indicators and the continuous improvement.



## SPECIFIC QC TESTS

With its extensive knowhow, the CARLO ERBA Reagents Quality Control team can not only manage and deliver the required tests needed for your product but also follow the requirements of your own test laboratory.



## CERTIFICATIONS

In order to always be in tune with the needs of the Pharmaceutical Industry market and to provide the best service from a quality point of view, we have been inspected and received the **ANSM certification** for excipients for our 2 production sites.



## CUSTOM PACKAGING

Produced regularly or following a custom request, we pack your chemicals in the container of your choice and then test it in accordance with the requirements of our quality system.


- Bottle, barrels, bucket with total opening, reusable container in polypropylene, glass, stainless steel, metal, from a few millilitres or milligrams to several tons.
- Inorganic or organic products, solid or liquid.
- Custom labelled, preweighing, pre-proportioning.



## SHUTTLE SERVICE

Stainless steel returnable containers to optimize solvent quality and the management of packaging waste

- Improved safety



- Environmental impact

- Quality preserved

- Efficient logistic

## DELIVERIES IN LARGE VOLUMES

The solvents you used in large quantities can be delivered in IBC and tankers from 500 to 30 000l. The product type and usage determine the type of IBC of the material.

Such deliveries will include:

- Certificate of analysis of the container
- Representative sample delivered with the IBC under request
- Computerized tracing of the maintenance and transport
- Full traceability of the product



# PACKAGING



## CARLO ERBA Reagents

Your partner in choice  
for laboratory and industry

With dedication, competence and experience, CARLO ERBA Reagents strives to be the right "partner in choice" offering products of certified quality to professionals involved with Research, Healthcare and Industry.



**CARLO ERBA**

REAGENTS

CARLO ERBA Reagents S.A.S  
CHAUSSÉE DU VEXIN  
27106 VAL DE REUIL - FRANCE  
TEL. +33 2 32 09 20 00  
FAX + 33 2 32 09 20 20

[www.carloerbareagents.com](http://www.carloerbareagents.com)  
[servicclient@cer.dgroup.it](mailto:servicclient@cer.dgroup.it)

